



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
REGION 10 LABORATORY
7411 Beach Dr. East
PortOrchard, Washington 98366

MAR 05 2015

MEMORANDUM

TO: Julie Wroble, Toxicologist
Risk Evaluation Unit
Office of Environmental Assessment
U.S. EPA Region 10

FROM: Jed Januch, Environmental Protection Specialist
Environmental Services Unit
Office of Environmental Assessment
U.S. EPA Region 10

SUBJECT: Case narrative for Technical Support – Fluidized Bed Asbestos Segregator
Sumas Mountain Asbestos Methods Study
Project Code: SFP-078A

A handwritten signature in black ink, appearing to read "Jed Januch".

Introduction

This memorandum documents sampling of soil materials conducted by staff from the U.S. EPA Region 10 Office of Environmental Assessment (OEA) using a fluidized bed asbestos segregator (FBAS). Forty-five specimens of soil were submitted for sampling in the FBAS. The sampling was done between January 26, 2015, and February 5, 2015, by Jed Januch, U.S. EPA Region 10.

Soil Preparation

The forty-five specimens of soil were dried in a laboratory oven at 60° centigrade (C) for 24 hours. After drying, the samples were sieved through a Number 20, USA Standard Testing Sieve with 850 micrometer (μm) mesh. The fraction greater than 850 μm was examined with the aid of a Wild M5 stereomicroscope to assess whether suspected asbestos structures were present. There were suspected asbestos structures detected in the fractions greater than 850 μm for some of the samples containing grey colored sand/sediment. The fractions less than 850 μm were tested in the FBAS.

Sampling Procedure

Sampling was done according to the U.S. EPA Region 10 standard operating procedure (SOP) OEAFIELDSOP-102 for Sampling, Sample Preparation and Operation of the Fluidized Bed Asbestos Segregator, Revision 1.0, dated August 19, 2011. To begin with, aliquots of each specimen (approximately 0.5 -1.0 grams) were combined with laboratory grade Ottawa 20/30 sand (approximately 19.0-19.5 grams). The sample/sand mixture (bed) was placed inside a clean glass vessel and air was drawn through the specimens for a sampling duration of 3 minutes. When the air flow is increased to where the pressure drop through the bed equals the weight of the bed, the solid material begins to circulate and act as a fluid. Small particles elutriate from the bed and are drawn with the air flow through the top of the glass vessel where a portion of the

air/particulate mixture is split off at an average flow rate of 237 cubic centimeters per minute (cc/m) and drawn through a 0.8 μ m pore size mixed cellulose ester (MCE) filter. The particulates captured on the filter will be analyzed by transmission electron microscope (TEM) using the counting rules of the ISO10312 method. A summary of the sampling information is included in Table 1.

Table 1 – Sample information summary

Fluidized Bed Asbestos Segregator - Sampling Record
 Project: Sumas Mtn. Asbestos Methods Study
 Project Code: SFP-078A

Date	Time	Sample ID	Grams Sample	Grams Sand	Comments
1/26/2015	1600	14394100	0.555	19.50	Light brown soil - PMC/SEM/EDS duplicate
1/26/2015	1619	14394101	0.501	19.50	Light brown soil
1/26/2015	1627	14394102	0.500	19.50	Dark brown soil - vegetation
1/26/2015	1715	14394103	0.509	19.50	Dark brown soil - vegetation
1/26/2015	1730	14394104	0.510	19.50	Dark brown soil - vegetation
1/28/2015	1545	14394105	0.495	19.50	Dark brown soil - vegetation
1/28/2015	1610	14394106	0.500	19.50	Dark brown soil - vegetation
1/28/2015	1630	14394107	0.508	19.50	Light brown soil
1/28/2015	1647	14394108	0.499	19.50	Light brown soil
1/28/2015	1655	14394109	0.501	19.50	Light brown soil
1/29/2015	1102	14394110	0.506	19.50	Light brown soil - vegetation
1/29/2015	1153	14394111	0.514	19.50	Light brown soil - vegetation
1/29/2015	1400	14394112	0.511	19.50	Light brown soil - vegetation
1/29/2015	1430	14394113	0.520	19.50	Light brown soil - vegetation
1/29/2015	1530	14394114	0.504	19.50	Light brown soil - vegetation
1/29/2015	1549	14394115	1.001	19.00	Light brown soil - vegetation
1/29/2015	1630	14394116	1.049	19.00	Light brown soil - vegetation
1/29/2015	1649	14394117	1.003	19.00	Light brown soil
1/29/2015	1700	14394118	1.006	19.00	Light brown soil - vegetation
1/29/2015	1715	14394119	1.005	19.00	Light brown soil
2/2/2015	1330	14394120	1.005	19.00	Light brown soil
2/2/2015	1500	14394121	1.009	19.00	Light brown soil
2/2/2015	1647	14394122	1.008	19.00	Light brown soil
2/2/2015	1714	14394123	1.001	19.00	Light brown soil
2/2/2015	1724	14394124	1.009	19.00	Light brown soil
2/3/2015	0916	14394125	1.008	19.00	Light brown soil
2/3/2015	0927	14394126	1.004	19.00	Light brown soil
2/3/2015	0942	14394127	1.003	19.00	Light brown soil - vegetation
2/3/2015	1348	14394128	1.000	19.00	Light brown soil
2/3/2015	1425	14394129	1.000	19.00	Light brown/grey soil
2/3/2015	1536	14394130	1.004	19.00	Light brown/grey soil
2/3/2015	1629	14394131	1.008	19.00	grey sand/sediment - PCM/SEM/EDS duplicate

2/3/2015	1643	14394132	0.998	19.00	grey sand/sediment
2/3/2015	1658	14394133	1.003	19.00	grey sand/sediment
2/4/2015	1336	14394134	1.009	19.00	grey sand/sediment
2/4/2015	1447	14394135	1.003	19.00	grey sand/sediment
2/4/2015	1502	14394136	0.999	19.00	grey sand/sediment
2/4/2015	1617	14394137	1.004	19.00	grey sand/sediment
2/4/2015	1630	14394138	1.000	19.00	grey sand/sediment
2/4/2015	1434	14394139	1.003	19.00	grey sand/sediment
2/5/2015	1529	14394140	1.001	19.00	grey sand/sediment
2/5/2015	1549	14394141	1.005	19.00	grey sand/sediment
2/5/2015	1559	14394142	1.002	19.00	grey sand/sediment
2/5/2015	1610	14394143	0.999	19.00	grey sand/sediment
2/5/2015	1723	14394143	1.000	19.00	grey sand/sediment - QA Replicate
2/5/2015	1627	14394144	1.004	19.00	grey sand/sediment

Quality Assurance and Quality Control

Sand blanks and process blanks were collected during the course of the project. Analysis of the sand blanks and the process blanks will help determine if cross contamination is occurring, and if so, the degree to which it occurs. In addition, two filter lot blanks from SKC Lot Number 14723-7DE-036 were submitted for analysis. A copy of the quality certification for the MCE filters is appended to this report. Table 2 provides a summary of the quality control (QC) samples prepared for this project.

Table 2 – QC sample summary

Fluidized Bed Asbestos Segregator - Quality Control Sample Summary

Project: Sumas Mtn. Asbestos Methods Study

Project Code: SFP-078A

Date	Sample ID	QC Sample Type	Q _r	Comments
12/9/2014	14394189	Lot Blank	n/a	SKC Lot No. 14723-7DE-036
12/9/2014	14394190	Lot Blank	n/a	SKC Lot No. 14723-7DE-036
12/9/2014	14394191	Sand Blank	0.0125	ResTek Ottawa Sand Lot No. 904101-AE
1/29/2015	14394192	Process Blank	n/a	Sample open face in HEPA workstation
2/2/2015	14394193	Sand Blank	0.0122	ResTek Ottawa Sand Lot No. 904101-AE
2/3/2015	14394194	Process Blank	n/a	Sample open face in HEPA workstation
2/5/2015	14394195	Sand Blank	0.0122	ResTek Ottawa Sand Lot No. 904101-AE

The flow rate of air passing through the FBAS was calibrated with a Gilian Primary Flow standard, serial number 003507, equipped with a high-flow cell, serial number 0208000-H for the total flow rate in liters per minute, and a low-flow cell, serial number 002598-S for filter calibration in cc/m. The calibrator base was serviced and the flow cells were calibrated by the manufacturer December 9-10, 2014. Copies of the calibration certificates for the flow cells are appended to this report.

The specimens of soil used for this project were weighed inside a HEPA hood with a Mettler PM400 balance, serial number M73584. The balance was checked daily with the aid of a Troemner ASTM Class 1 metric weight set, serial number 41766. The balance was accurate within +/- 0.001 grams. The sand was weighed outside the HEPA hood with a Mettler Toledo B502-S balance serial number 1119200622. It was also checked with the Troemner metric weight set daily and found to be accurate within +/- 0.01 grams.

The degree of filter loading on two samples was screened with the aid of an Olympus CX21 phase contrast microscope (PCM) at 400x magnification. Images of random locations on the filters were collected using a Carl Zeiss Axioskop 40 microscope (set up for PCM) equipped with a Carl Zeiss AxioCam. In addition, two replicate samples were prepared in polycarbonate filters for analysis by scanning electron microscope (SEM) and energy dispersive spectroscopy (EDS). Photomicrographs by PCM and SEM as well as EDS spectra are appended to this report.

Sample Shipping and Chain of Custody

Samples were hand delivered by EPA staff on February 9, 2015, to Lab/Cor, Inc., located at 7619 6th Avenue NW, Seattle, WA 98117 and were received for Lab/Cor by Kate March. The original chain of custody forms are appended to this report.

Appendices

1. Quality Certificate – MCE filters
2. Flow Cell Calibration Certificates
3. Photomicrographs
4. Chain of Custody Forms

Appendix 1

CERTIFICATION

SKC Omega Specialty Division is committed to manufacturing the highest quality cassettes for asbestos sampling. This commitment is backed by the BestChek™ Seal of Approval which assures a higher standard in cassette reliability. Our BestChek™ approved cassettes are assembled in a full HEPA filtered clean environment. Our cassettes comply with both the NIOSH 7400 method which allows up to 5 fibers/100 fields and OSHA ID160 regulations which allow 4 fibers/100 fields. Always look for the BestChek™ Seal of Approval when purchasing air monitoring cassettes.

Cassette Category: 25mm, 3 Piece Conductive Cassette

Product Code: 225-326

Lot Number: 14723-7DE-036

Compliance: NIOSH 7400, OSHA ID-160

Average Background Count: 1.75 fibers/100 Fields

PCM Analysis: Passed

Clearing, Acetone Vapor: Passed

Clearing, DP/DO (Monitoring Solution): Passed

Plastic Conductivity by ASTM D-257: 2.3×10 ohms-cm

Percent Carbon by ASTM D-1603: 19.05

Cowl Conductivity by Resistance Meter: 10-5008 avg/ (x1000)

Filter Material: MCE

Pore Size: 0.8 μ m


Pad Material: Cellulosic

Effective Filter Area: 385 mm²

Cowl Length: 2 Inch

Analyst: PH

Please keep for your records



1-800-752-8472

Appendix 2

Sensidyne, LP

CALIBRATION CERTIFICATE

Cell S/N: 0208000-H

Date: December 10, 2014

This is to certify that the above referenced Gilibrator Flow Cell was calibrated using film flowmeter MCH-101-A, which has been calibrated by instruments directly traceable to the National Institute of Standards and Technology. NIST Report 8361604.

Results:

REFERENCE MCH-101-A l/min	S/N 0208000-H l/min	RELATIVE DIFF. l/min	PERCENT DIFF.
5.036	5.032	-0.004	-0.08
5.031	5.031	0.000	0.00
5.038	5.040	0.002	0.04
5.036	5.037	0.001	0.02
5.040	5.034	-0.006	-0.12
5.039	5.031	-0.008	-0.16
5.038	5.038	0.000	0.00
5.036	5.032	-0.004	-0.08
5.039	5.036	-0.003	-0.06
5.037	5.029	-0.008	-0.16
<hr/>			
MAX		-0.008	-0.160
MEAN	5.037	5.034	

Calibrated by

Agnes Barfield

Date: December 10, 2014

CODE 000

Sensidyne, LP

CALIBRATION CERTIFICATE

Cell S/N: 002598-S

Date: December 10, 2014

This is to certify that the above referenced Gilibrator Flow Cell was calibrated using film flowmeter MCS-102-A, which has been calibrated by instruments directly traceable to the National Institute of Standards and Technology. NIST Report 8361604.
Results:

REFERENCE MCS-102-A cc/min	S/N 002598-S cc/min	RELATIVE DIFF. cc/min	PERCENT DIFF.
2017	2018	1	0.05
2016	2018	2	0.10
2018	2018	0	0.00
2016	2020	4	0.20
2021	2017	-4	-0.20
2020	2018	-2	-0.10
2018	2018	0	0.00
2020	2018	-2	-0.10
2022	2018	-4	-0.20
2020	2020	0	0.00

MAX

4.0

0.2

MEAN 2018.8

2018.3

Calibrated by

Agnes Beal

Date: December 10, 2014

CODE 100

Appendix 3

Sumas Mountain Asbestos Methods FBAS Test Samples Project Code: SFP-078A

Analyst

Jed Januch

Equipment

Stereomicroscope – Wild M5 with Nikon CoolPix camera

PLM – Carl Zeiss Axioskop 40 with AxioCam MRc imaging system

SEM – JEOL JSM6510 LV

EDS – EDAX Genesis System

Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A

Sumas Mtn Asbestos Methods Study
Project Code: SFP-078A
Sample: 14394100
FBAS 1 gram sample
400x, PCM
Jed Januch, 2/25/15

100 μ m

Sumas Mtn Asbestos Methods Study
Project Code: SFP-078A
Sample: 14394100
FBAS 1 gram sample
400x, PCM
Jed Januch, 2/25/15

100 μ m

Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A

Sumas Min Asbestos Methods Study
Project Code: SFP-078A
Sample: 14394100
FBAS 1 gram sample
400x, PCM
Jed Januch, 02/24/15

20 μ m

Sumas Min Asbestos Methods Study
Project Code: SFP-078A
Sample 14394131
FBAS, 1 gram sample
400x, PCM
Jed Januch, 3/2/2015

20 μ m

Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A


Sumas Mtn Asbestos Methods Study
Project Code: SFP-078A
Sample 14394131
FBAS, 1 gram sample
400x, PCM
Jed Januch, 3/2/2015

20 μ m



Sumas Mtn Asbestos Methods Study
Project Code: SFP-078A
Sample 14394131
FBAS, 1 gram sample
400x, PCM
Jed Januch, 3/2/2015

20 μ m



Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A





Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A





Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A





Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A





Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A



Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A



Sumas Mountain Asbestos Methods
FBAS Test Samples
Project Code: SFP-078A



Appendix 4

EPA Region 10 - FBAS CHAIN OF CUSTODY RECORD
OEAFIELD SOP-102 (Revision 1.0)

COC No.: _____
[FBAS-YYMMDD###]

COC entered by (signature):

COC Page 1 of 1

Project Manager: Julie Winkle

Project Identifier: Synas Mtn. Asbestos - FBAS

Method of Shipment: *EPA - deliver*

Shipment Tracking No.:

Ship Date:

Shipment Destination: Lab/Cov, Inc.
Seattle, WA

Shipment Origin: EPA Region 10 Laboratory
Port Orchard, WA

[illegible]

⁺⁺ Prepare replicate filters using fluidized bed asbestos segregator in accordance with SOP.

** Specific recording and stopping rules are specified in a site/project specific QAPP. As needed, grids may be prepared indirectly in accordance with ISO 13794:1999(E).

FBAS CHAIN OF CUSTODY RECORD

COC No.: FBAS-150209001

[FBAS-YYMMDD###]

COC entered by (signature): Julie JohnCOC Page 1 of 3Project Manager: Julie JohnProject Identifier: Sumas Mtn Asbestos Methods - FBASMethod of Shipment: Hand DeliveredShipment Tracking No.: N/AShip Date: 2/9/15Shipment Origin: EPA Region 10 Laboratory
Port Orchard, WAShipment Destination: LABCOR, Inc.
Seattle, WA

Sample Information		Sample Preparation Request			Filter Analysis Request				FBAS Information			Sample Remarks
Sample ID (PE Standard ID or Replicate Filter ID)	Medium [S for soil; A for air filter]	OEAFIELD SOP- 102 (Rev 1)** with cyclone	without cyclone		TEM, ISO 10312:1995(E)**				Blank [LT for lot; PB for prep; SB for sand]	FBAS Flow Ratio (Q _R)	FBAS Soil Mass (grams, dry weight)	
14394100	A	X	X	X						0.0122	0.555 g	1/27/15
14394101	A	X	X	X						0.0122	0.501 g	1/27/15
14394102	A	X	X	X						0.0122	0.500 g	1/27/15
14394103	A	X	X	X						0.0122	0.509 g	1/27/15
14394104	A	X	X	X						0.0122	0.510 g	1/27/15
14394105	A	X	X	X						0.0122	0.495 g	1/28/15
14394106	A	X	X	X						0.0122	0.500 g	1/28/15
14394107	A	X	X	X						0.0122	0.508 g	1/28/15
14394108	A	X	X	X						0.0122	0.499 g	1/28/15
14394109	A	X	X	X						0.0122	0.561 g	1/28/15
14394110	A	X	X	X						0.0122	0.506 g	1/29/15
14394111	A	X	X	X						0.0122	0.514 g	1/29/15
14394112	A	X	X	X						0.0122	0.511 g	1/29/15
14394113	A	X	X	X						0.0122	0.520 g	1/29/15
14394114	A	X	X	X						0.0122	0.504 g	1/29/15
14394115	A	X	X	X						0.0122	1.001 g	1/29/15
14394116	A	X	X	X						0.0122	1.049 g	1/29/15
14394117	A	X	X	X						0.0122	1.003 g	1/29/15
14394118	A	X	X	X						0.0122	1.006 g	1/29/15
14394119	A	X	X	X						0.0122	1.005 g	1/29/15
14394120	A	X	X	X						0.0122	1.005 g	2/2/15

RELINQUISHED BY:

NAME, AFFILIATION

DATE, TIME

RELINQUISHED TO:

NAME, AFFILIATION

DATE, TIME

Julie John US: EPA

2/9/15 1330

Julie John YC

2/9/15 1330

COMMENTS/CONDITION OF SAMPLES:

** Prepare replicate filters using fluidized bed asbestos segregator in accordance with SOP.

** Specific recording and stopping rules are specified in the study QAPP. As needed, grids may be prepared indirectly in accordance with ISO 13794:1999(E).

FBAS CHAIN OF CUSTODY RECORD

COC No.: FBAS-150209002

(FBAS-YYMMDD###)

COC entered by (signature): [Signature]COC Page 2 of 3Project Manager: Julie WohleProject Identifier: Seas Man Asbestos Methods - FBASMethod of Shipment: Hand ReliefShipment Tracking No.: N/AShip Date: 2/9/15Shipment Origin: EPA Region 10 LaboratoryShipment Destination: WAB/CR, Inc.Port Orchard, WASeattle, WA

Sample Information		Sample Preparation Request			Filter Analysis Request			FBAS Information		Sample Remarks
Sample ID (PE Standard ID or Replicate Filter ID)	Medium [S for soil; A for air filter]	ORAFIELD SOP- 102 (Rev 1)** with cyclone	without cyclone	TEM, ISO 10312:1999(E)**	Blank [LT for lot; PB for prep; SB for sand]	FBAS Flow Ratio (Q _a)	FBAS Soil Mass (grams, dry weight)			
14394121	A	X	X	X		0.022	1.009 g	2/2/15		
14394122	A	X	X	X		0.022	1.008 g	2/2/15		
14394123	A	X	X	X		0.022	1.001 g	2/2/15		
14394124	A	X	X	X		0.022	1.009 g	2/2/15		
14394125	A	X	X	X		0.022	1.008 g	2/3/15		
14394126	A	X	X	X		0.022	1.004 g	2/3/15		
14394127	A	X	X	X		0.022	1.003 g	2/3/15		
14394128	A	X	X	X		0.022	1.000 g	2/3/15		
14394129	A	X	X	X		0.022	1.000 g	2/3/15		
14394130	A	X	X	X		0.022	1.004 g	2/3/15		
14394131	A	X	X	X		0.022	1.008 g	2/3/15		
14394132	A	X	X	X		0.022	0.998 g	2/3/15		
14394133	A	X	X	X		0.022	1.003 g	2/3/15		
14394134	A	X	X	X		0.022	1.009 g	2/4/15		
14394135	A	X	X	X		0.022	1.003 g	2/4/15		
14394136	A	X	X	X		0.022	0.999 g	2/4/15		
14394137	A	X	X	X		0.022	1.004 g	2/4/15		
14394138	A	X	X	X		0.022	1.000 g	2/4/15		
14394139	A	X	X	X		0.022	1.003 g	2/4/15		
14394140	A	X	X	X		0.022	1.001 g	2/5/15		
14394141	A	X	X	X		0.022	1.005 g	2/5/15		

RELINQUISHED BY:		RELINQUISHED TO:	
NAME, AFFILIATION	DATE, TIME	NAME, AFFILIATION	DATE, TIME
<u>[Signature]</u> J.S.EPA	2/9/15 1330	<u>Wohle</u> L/C	2/9/15 1330

COMMENTS/CONDITION OF SAMPLES:

** Prepare replicate filters using fluidized bed asbestos segregator in accordance with SOP.

** Specific recording and stopping rules are specified in the study QAPP. As needed, grids may be prepared indirectly in accordance with ISO 13794:1999(E).

FBAS CHAIN OF CUSTODY RECORD

COC No.: FBAS-150209003

COC entered by (signature):

[FBAS-YYMMDD###]
COC Page 3 of 3

Project Manager: Julie Washburn

Project Identifier: Sinus with Asbestos Methods - FBAS

Method of Shipment: Hand Delivered

Shipment Tracking No.: N/A

Ship Date: 2/9/15

Shipment Origin: EPA Region 10 Laboratory

Shipment Destination: LAB/COR, Inc.

Port Orchard, WA

Seattle WA

[illegible]

^{**} Prepare replicate filters using fluidized bed asbestos segregator in accordance with SOP.

** Specific recording and stopping rules are specified in the study QAPP. As needed, grids may be prepared indirectly in accordance with ISO 13794:1999(E).